



# PRODUCT CHOICE





MediWipes nonwoven wipes.....	6-9
MediSept rapid sensitiv.....	10/11
MediSept uni.....	12/13
Sachets.....	14/15
Quick & Easy rapid sensitiv.....	16/17
Quick & Easy uni.....	18/19
MediWipes dispenser systems.....	20/21

SmartPoint.....	22/23
Smart Bag.....	24/25

FloorWipes.....	26/27
-----------------	-------

PowerWipes.....	28/29
clean&care Poly-Pro Wipes.....	30/31
clean&care Vis-Co Wipes.....	32/33

Wall brackets.....	34/35
--------------------	-------

MediWash.....	36/37
---------------	-------



The top principle of our company's philosophy is to meet the individual requirements of our customers. In order to do this we provide not only high quality, but also user friendly and efficient products at the highest quality levels. Since its foundation in 1999, we have developed into a highly qualified producer of fleece and achieved a leading position in the market.

Well thought out warehouse and shipping logistics guarantee reliable and quick delivery. You will find our products both nationally and internationally, mostly in private label packaging. Customer-oriented developed and assembled, they meet all the requirements of functionality, design and efficiency.

Absolute reliability in quality, adherence to deadlines and sustainability are our commitment and for us, a matter of course. This means concentrated competence which we demonstrate on a daily basis.

Our products are characterised by „MADE IN GERMANY“ quality at fair prices. We enjoy the great confidence that our customers have in us and hold true to our responsibility in maintaining this.

We see ourselves not only as a reliable manufacturer and supplier, but also as a provider of inspiration and co-creation. Innovativeness, flexibility, quality, customer focus and punctuality are traits appreciated by our clients. We rely on partnerships and team work.



# MediWipes

## Disposable wipes for the hygienic disinfection of surfaces

### Non-woven wipes roll

The MediWipes consist of a 100% Polyester Fleece and ensures an optimal substance release during application. Also, because of the cross-linking, they are extremely tear-resistant and hard-wearing.

The MediWipes disinfection process is ready for use. Each non-woven wipes roll is labelled with a batch label which documents the batch, production- and expiry date of the disinfectant solution.

The single wipe removal from the re-sealable dispenser system means the MediWipes are used very economically.

### Application



1. Put on gloves.



2. Remove the non-woven wipes roll from the packaging and insert into the purified, disinfected and dry wet wipes dispenser.



3. Wet the non-woven wipes roll with the recommended amount.



4. Lead the non-woven wipes through the opening from beneath.



5. Fill in the attached batch label and stick on.

## Variants

### MediWipes Type I

(Article number see "Delivery Forms")

The non-woven wipes rolls are individually wrapped in foil and shrink-wrapped



+

### MediWipes Type II

(Article number see "Delivery Forms" + suffix "SA")

The non-woven wipes rolls are individually wrapped in foil and shrink-wrapped

+ WHR cleaning sachet  
+ WHR cleaning-sachet



+

### MediWipes Type III

(Article number see "Delivery Forms" + suffix "AI")

The non-woven wipes rolls are individually wrapped in foil and shrink-wrapped

+ WHR cleaning sachet  
+ WHR cleaning-sachet  
+ Foil bag and attachment ring  
+ Timestrip time indicator



# MediWipes

## Disposable wipes for the hygienic disinfection of surfaces

### Delivery forms

Wipe quality	Material	Wipe size	Number of wipes	Article no.	VPE	Dispenser system	
premium	PET <sub>smooth</sub>	approx. 30 x 30 cm	30	201030302601	6 rolls	can „mini“	
			30	20103030	6 rolls + 1 can	can „mini“	
			60	201060302601	8 rolls	bucket „midi“	
			70	201070302601	8 rolls	bucket „midi“	
			80	201080302601	8 rolls	bucket „midi“	
		approx. 30 x 32 cm	90	201090302601	6 rolls	bucket „maxi“	
			50	201050322601	4 rolls	bucket „midi“	
			100	201100322601	6 rolls	bucket „maxi“	
			approx. 20 x 36 cm	40	201040362602	5 rolls	can „XXL“
				90	201090362602	6 rolls	bucket „XXL“

Wipe quality	Material	Wipe size	Number of wipes	Article no.	VPE	Dispenser system
eco	PET <sub>smooth</sub>	approx. 20 x 36 cm	40	201040362402	5 rolls	can „XXL“
			90	201090362402	6 rolls	bucket „XXL“

Wipe quality	Material	Wipe size	Number of wipes	Article no.	VPE	Dispenser system
basic	PET <sub>perforated</sub>	approx. 30 x 32 cm	50	201050321501	10 rolls	bucket „midi“
			90	201090321501	8 rolls	bucket „maxi“

Wipe quality	Material	Wipe size	Number of wipes	Article no.	VPE	Dispenser system
plus	PET/VIS <sub>smooth</sub>	approx. 30 x 30 cm	70	203070302501	8 rolls	bucket „midi“
			90	203090302501	8 rolls	bucket „maxi“

\* You can find the suitable dispenser systems on pages 20 & 21.

## Ready to use quick disinfection

### Preparation

MediSept rapid sensitiv is a ready-to-use disinfectant, based on alcohol and on quaternary ammonium compounds. Serving for the comprehensive and quick disinfection of surfaces and objects. The cleaning agent is free of aldehydes, phenol, chlorine and alkali.

### Application

MediSept rapid sensitiv can be put on a disposable wipe. Alternatively, a pre-moistened role is also available. A corresponding data sheet is available. Thereupon, the area to be disinfected is completely and evenly wetted with the ready-to-use solution. Subsequently, the cleaning agent must work and the disinfected surfaces must dry. The working times are listed under the same named point (see below). There are almost no surface residues formed.

### Range of effect

bactericidal  
yeastcid  
tuberculocide  
limited virus-deactivating

### Application time

Noro viruses	5 min.
HCV	30 sec.
HBV, HIV	30 sec.
Vaccinia	30 sec.
TB bacteria	2 min.
Bacteria incl. MRSA	1 min.
Yeast fungus	1 min.

## Product status / listings

The MediSept rapid sensitiv impregnation liquid is included in the current list of disinfectants of the German Association for Applied Hygiene (VAH). Registered as a biocidal product at the BAuA (Federal Institute for Occupational Safety and Health) in accordance with the German Biocide Act: N-71152

## Delivery forms

Name	Containers	Article-No.	VPE
MediSept rapid sensitiv	1 litre	MS10-01	6 bottles
	10 litre	MS10-10	1 canister



BAuA-Nr. N-71233



## Alcohol free quick disinfection

### Preparation

MediSept uni is a alcohol free ready-to-use disinfectant. It serves for comprehensive quick disinfection of surfaces and objects, including those that are alcohol-sensitive. The cleaning agent is free of aldehydes, phenol, chlorine and alkali.

### Application

MediSept uni is put on a disposable wipe. Alternatively, a pre-moistened role is also available (quick & easy uni, Art.-No. MSW20500u, pages 16 & 17. Thereupon, the area to be disinfected is completely and evenly wetted with the ready-to-use solution. Subsequently, the cleaning agent must work and the disinfected surfaces must dry. The working times are listed under the same named point (see below). Following the working time, if necessary, wipe again with a disposable wipe. The alcohol free quick disinfection is odourless.

### Range of effect

bactericidal incl. MRSA  
yeastocidal  
virus-inactivating (HBV / HIV, HCV, rota, papova- /  
polyoma viruses)

### Application time

bactericidal incl. MRSA,  
yeastocidal according to the German Association for Applied Hygiene (VAH), under high  
organic contamination) 1 min.  
effective against all enveloped  
viruses including vaccinia, BVDV, HBV,  
HIV, HCV, influenza 30 sec.  
Rota viruses 30 sec.  
papova- / polyoma viruses 2 min.

## Product status / listings

Included in the current list of disinfectants of the German Association for Applied Hygiene (VAH) and the Austrian Society for Hygiene and Microbiology and Preventative Medicine (ÖGHMP).

Registered as a biocidal product at the BAuA (Federal Institute for Occupational Safety and Health) with the German Biocide Act: N-58067

## Delivery forms

Name	Containers	Article no.	VPE
MediSept uni	1 litre	MS20-01	12 bottles
	10 litre	MS20-10	1 canister



BAuA-Nr. N-58067

# Sachet

## Ready to use disinfection wipe

### alcohol-containing

#### Preparation

The ready-to-use wipes serve for comprehensive quick disinfection of alcohol-resistant surfaces and objects. They are saturated with MediSept rapid, a ready-to-use alcohol-based disinfectant with a melon fragrance. The cleaning agent is free of aldehydes, phenol, chlorine and alkali.

Best possible agent release is guaranteed through the high quality of the non-woven fabric. The wipes are packed separate.

#### Application

The area to be disinfected is completely and evenly wetted with the ready-to-use solution. Subsequently, the cleaning agent must work and the disinfected surfaces must dry. The working times are listed under the same named point (see right). No surface residues are formed. This product is not to be used with acrylic glass (Plexiglas).

#### Range of effect

The data refers to the impregnation liquid,

MediSept rapid:

bactericidal

fungicide

tuberculocide

limited virus-deactivating

#### Application time

The data refers to the impregnation liquid,

MediSept rapid:

Noro viruses (low impact) 30 sec.

Rota viruses 15 sec.

HCV, BVDV 15 sec.

HBV, HIV 5 min.

Vaccinia 15 sec.

Enveloped viruses 15 sec.

TB

(with mechanics, high impact) 30 sec.

bactericidal 30 sec.

fungicide 30 sec.

## Delivery forms

Name	Box	Article no.	VPE
Sachet	10 wipes	SD2018	18 boxes



# QUICK & EASY

*rapid sensitiv*

## Ready for use disinfectant wipes

### alcohol-containing

#### Preparation

The ready-to-use wipes serve for comprehensive quick disinfection of alcohol-resistant surfaces and objects. They are impregnated with MediSept rapid sensitiv, a ready-to-use alcohol and quaternary ammonium compounds based disinfectant. The cleaning agent is free of aldehydes, phenol, chlorine and alkali.

Best possible agent release is guaranteed through the high quality of the non-woven fabric.

#### Application

The non-woven wipes roll is placed into the dispenser. A disinfectant wipe is then led from the center of the roll through the opening of the lid. The lid is then closed again.

The area to be disinfected is completely and evenly wetted with the ready-to-use solution. Subsequently, the cleaning agent must work and the disinfected surfaces must dry. The working times are listed under the same named point (see right). There are almost no surface residues formed.

#### Range of effect

The data refers to the impregnation liquid,

MediSept rapid sensitiv:

bactericidal

yeastcid

tuberculocide

limited virus-deactivating

#### Application time

The data refers to the impregnation liquid,

MediSept rapid sensitiv:

Noro viruses 5 min.

HCV 30 sec.

HBV, HIV 30 sec.

Vaccinia 30 sec.

TB bacteria 2 min.

Bacteria incl. MRSA 1 min.

Yeast fungus 1 min.

## Product status / listings

The MediSept rapid sensitiv impregnation liquid is included in the current list of disinfectants of the German Association for Applied Hygiene (VAH). Registered as a biocidal product at the BAuA (Federal Institute for Occupational Safety and Health) in accordance with the German Biocide Act: N-71152

## Delivery forms

Name	Wipe size	Number of wipes	Article no.	VPE
quick & easy rapid sensitiv	approx. 30 x 30 cm	30 leaves	MSW2030	6 rolls
		70 leaves	MSW2070	3 rolls
	approx. 30 x 15 cm	60 leaves	6000-2003	6 rolls



BAuA-Nr. N-71152

\* You can find the suitable dispenser systems on pages 20 & 21.

## Ready for use disinfectant wipes

### alcohol-free

#### Preparation

The ready-to-use wipes serve for comprehensive quick disinfection of surfaces and objects. They are impregnated with MediSept uni, a ready-to-use disinfectant based on quaternary ammonium compounds. The cleaning agent is free of aldehydes, phenol, chlorine and alkali. Best possible agent release is guaranteed through the high quality of the non-woven fabric.

#### Application

The non-woven wipes roll is placed into the dispenser. A disinfectant wipe is then led from the center of the roll through the opening of the lid. The lid is then closed again.

The area to be disinfected is completely and evenly wetted with the ready-to-use solution. Subsequently, the cleaning agent must work and the disinfected surfaces must dry. The working times are listed under the same named point. Following the working time, if necessary, wipe again with a disposable wipe. The quick disinfection and cleaning is odourless.

#### Range of effect

The data refers to the MediSept uni impregnation liquid:  
bactericidal incl. MRSA  
yeastocidal  
virus-inactivating (HBV / HIV, HCV, rota, papova- / polyoma viruses)

#### Application time

The data refers to the MediSept uni impregnation liquid:  
bactericidal incl. MRSA,  
yeastocidal according to the German Association for Applied Hygiene (VAH), under high organic contamination) 1 min.  
effective against all enveloped viruses including vaccinia, BVDV, HBV, HIV, HCV, influenza 30 sec.  
Rota viruses 30 sec.  
Papova- / polyoma viruses 2 min.

## Product status / listings

The MediSept uni impregnation liquid is included in the current list of disinfectants of the German Association for Applied Hygiene (VAH) and the Austrian Society for Hygiene and Microbiology and Preventative Medicine (ÖGHMP). Registered as a biocidal product at the BAuA (Federal Institute for Occupational Safety and Health) with the German Biocide Act: N-58068

## Delivery forms

Name	Wipe size	Number of wipes	Article no.	VPE
quick & easy uni	approx. 30 x 30 cm	30 leaves	MSW2030u	6 rolls
		70 leaves	MSW2070u	3 rolls
	approx. 20 x 20 cm	500 leaves	MSW20500u	2 rolls



BAuA-Nr. N-68231

\* You can find the suitable dispenser systems on pages 20 & 21.

# MediWipes

## Dispenser systems for wipes for the hygienic disinfection of surfaces

### Dispenser systems

The MediWipes dispenser consist of polyethylene and polypropylene. The bucket / can is therefore resistant to all conventional disinfectants and cleaning agents.

The wipes can be easily removed through the opening. The star-shaped indentation in the pre-moistened wipes opening system prevents the wipes from drying out. The wipes are available in different colours. These are listed under „variants“. The transparent lid of the buckets allows for easy monitoring of the pre-moistened wipes.

The bucket / the can and the corresponding lid are dishwasher-safe up to min. 65°C.

### Application



1. Put on gloves.



2. Remove the non-woven wipes roll from the packaging and insert into the purified, disinfected and dry wet wipes dispenser.



3. Wet the non-woven wipes roll with the recommended amount.



4. Lead the non-woven wipes through the opening from beneath.



5. Fill in the attached batch label and stick on.

## Variants

The pre-moistened wipes opening system is available in different colours, to allow a distinction between application areas, used agents etc. The following colours are available:

## Delivery forms

Dispenser system	Large	Volumes	Article no.	VPE
dispenser	mini	approx. 1.7 Litre	1011700V	5 cans
	XXL	approx. 2.5 Litre	1012000V	5 cans
dispenser bucket	mini	approx. 2.6 Litre	10226175V	12 buckets
	midi	approx. 3.8 Litre	10238200V	8 buckets
	maxi	approx. 5.6 Litre	10256225V	6 buckets
	XXL	approx. 6.2 Litre	10262225V	5 buckets
	quattro	approx. 5 Litre	102503550	5 buckets

\* You can find the suitable wall brackets on pages 30 & 35.



# SMARTPOINT

## Free-standing dispenser system

### Dispenser System

The SmartPoint is a free-standing disinfectant wipes dispenser which can be placed anywhere. It is available in different colours and can be used for many different areas of application.

The used disinfectant wipes can be immediately discarded via the integrated waste management system.

### Application

In the SmartPoint fits the quick & easy uni roll (Art.-No. MSW20500U, pages 16 & 17). The first wipe must be fed through the opening from the inside. The SmartPoint is then ready for use.

## Delivery forms

Product	Article no.	Colour	
SmartPoint	SP105	green	
		yellow	
		orange	
		red	
		blue	
		black	
		granite	
		sandstone	
		white marble	
		SP106	stainless steel



# SMARTBAG

## Disposable wipes for the hygienic disinfection of surfaces

### Non-woven wipes roll


The SmartBag-System consists of a dry non-woven wipes roll and a stand-up bag with carrying handle. Each stand-up bag is equipped with a pre-moistened wipes opening system. This extraction can be twisted out, allowing for wetting of the wipes roll with any surface disinfectant. The non-woven wipes roll consists of a 100% Polyester Fleece and ensures an optimal substance release during application. Also, because of the cross-linking, it is extremely tear-resistant and hard-wearing.

The single wipe removal from the re-sealable system means the SmartBag is very economical.

### Application



1. Put on gloves.



2. Remove the tamper-proof seal.



3. Remove the opening system by turning counter-clockwise.



4. Add approx. 2.5 litres of ready-to-use disinfectant solution.



5. Lead the first wipe through the back of the opening system.



6. Insert the opening system again by turning in a clockwise direction.



7. Thoroughly wipe the surfaces or devices, while ensuring complete wetting and adhere to the working time of the disinfectant.

## Areas of application

The SmartBag system is specifically designed for areas with high risk of infection. The dispenser's one-way system minimises the risk of cross-contamination. Disposal costs are minimised.

## Delivery forms

Wipe quality	Material	Wipe size	Number of wipes	Article no.	VPE
on request					



# FloorWipes

## Disposable wipes for the hygienic floor disinfection

### Non-woven wipes roll

The FloorWipes consist of a 100% Polyester Fleece and ensures an optimal substance release during application. Also, because of the cross-linking, they are extremely tear-resistant and hard-wearing.

The FloorWipes wipe system is ready for use. Each non-woven wipes roll is labelled with a batch label which documents the batch, production- and expiry date of the disinfectant solution.

The single wipe removal from the re-sealable dispenser system means the FloorWipes are used very economically.

### Application



1. Put on gloves.



2. Remove the FloorWipes non-woven wipes roll out of the packaging and place in an appropriate, disinfected and dry dispenser system.



3. Depending on the roll size, sufficient ready-to-use disinfectant solution must be added.



4. Lead the non-woven wipes through the opening from beneath.



5. Fill in the attached batch label and stick on.



6. Remove the wipe cloth, lay out on the floor and attach using hook and loop connection to the mop board.



7. Thoroughly mop the floor, while ensuring complete wetting and adhere to the working time of the disinfectant.

## Delivery forms

Wipe quality	Material	Wipes size	Number of wipes	Article no.	VPE	Dispenser system
premium	PET <small>smooth</small>	approx. 20 x 44 cm	120	2011204420602	4 Rollen	Eimer „XXL“

## Accessories

Name	Containers	Article no.	VPE
Package Bracket 40	Mop holder 40 cm wide + aluminium telescopic handle	410301	1 holder + 1 handle

\* You can find the suitable wall brackets on pages 34 & 35.





## Saturated Cleaning Wipes

### Liquid

The ready-to-use wipes serve for the cleaning of dirty hands, machines and tools. They are impregnated with a pleasantly scented hand cleaner that contains skin-friendly cleansing agents.

Hand cleaning is made possible in areas where is no available soap and water.

The 3 product variants differ in the specifications of the non-woven material.

### PowerWipes

The non-woven consists of a synthetical fibre mixture. The **white** scarf is absorbent and collects the dirt. They can be used by heavy contamination of fats, grease, oils, ink, wax, tar, coal, graphite, grass stains and so on.

### PowerWipes Pro

The non-woven consists of 100% polypropylene. The **blue** scarf is very absorbent and can be used several times (squeezable). The PowerWipes Pro can be used for cleaning of machines and tools, because they are poor of lints and collect fats and oils in an optimally.

### PowerWipes Xtrem

The non-woven consists of 100% polypropylene. The **orange** scarf is abrasive and able to remove also heavy containments. The PowerWipes Xtrem scrub very dirty hands and machines.

Further they can be used by heavy contamination of fats, grease, oils, ink, wax, tar, coal, graphite, grass stains and so on.

## Application

Remove a single wipe from the opening of the dispenser system. Thoroughly clean hands / surfaces with the wipe.

## Delivery forms

Name	Option	Number of wipes	Wipe size	Article-No.	VPE
PowerWipes		75 leaves	approx. 20 x 36 cm	PW2075	6 cans
	Refill			PW2075R	6 cans
	Starter kit			PW2075S	1 can + 5 rolls
PowerWipes Pro		75 leaves	approx. 20 x 30 cm	PW2075P	6 cans
	Refill			PW2075PR	6 rolls
PowerWipes Xtrem		75 leaves	approx. 20 x 36 cm	PW2075X	6 cans
	Refill			PW2075XR	6 rolls



\* You can find the suitable wall brackets on pages 34 & 35.

## Poly-Pro Wipes

### Non-woven

The scarf can be used for cleaning of dirty machines and tools. It is wet strength, absorbent and almost lint-free. Further the Poly-Pro Wipes are solvent resistant and collect oil. The wipes are squeezable and can be used several times. The blue non-woven material consists of 100 % polypropylene and has a star-shaped structure.

### Application

Remove a single wipe. The dirty surface can be cleaned with the dry wipe or also with moistened wipe.

## Areas of application

The Poly-Pro Wipes have been proved of value in garages and industry companies. Also oily contaminations can be removed with them. They can be used for cleaning print rolls and for cleaning in colour- and varnish-using companies, because they are solvent resistant. The Poly-Pro wipes are also proved for canteen kitchens, because they are food save.

## Delivery forms

Name	Number of wipes	Article-No.	Palette
Clean & Care Poly-Pro Wipes	400 leaves	302400	48 rolls

A suitable floor stand is available.



## Vis-Co Wipes

### Non-woven

The scarf can be used for cleaning of dirty machines and tools. It can absorb big amounts of liquid and stays tear-resistant. Vis-Co Wipes are poor of lints, solvent resistant and food safe. The white non-woven material consists of 100% viscose and has a fine hole-structure.

### Application

Remove a single wipe. The dirty surface can be cleaned with the dry wipe or also with moistened wipe.

## Areas of application

The Vis-Co Wipes have been proved of value in garages and industry companies. They can be used for cleaning print rolls and for cleaning in colour- and varnish-using companies, because they are solvent resistant. Because of the finde hole-structure, polish kiquids and wax can be absorbed. The Vis-Co wipes are also proved for canteen kitchens, because they are food save.

## Delivery forms

Name	Wipes size	Number of wipes	Article no.	VPE	Pallete
Clean & Care Vis-Co Wipes	approx. 40 x 38 cm	400 leaves	305400	1 rolls	48 rolls
	approx. 28 x 38 cm	50 leaves	3030503716031	12 pack (z-folded)	48 VPE

A suitable floor stand is available.



# Wall-brackets

## Fixings for every application

### Wall bracket

The practical, modern design of the wall brackets are used for the safe keeping of all dispenser buckets and cans. The wall bracket is easy to install with little effort. Fixing materials are included. Due to the smooth surface, the wall mount is easy to disinfect.

### Application

The wall brackets are fixed with the supplied fixing material at the desired location. The matching buckets / cans are then set in the holder. Depending on the variant, the wall holder can still be locked in order to prevent unwanted removal of the dispenser systems by unauthorised persons.

### Variants

#### Eco-Series

The wall brackets are made of powder coated steel and are white. A suitable wall bracket is available for each bucket size. There are 2 different variants of which one is lockable ("eco" safe) and one is not ("eco" basic). The „eco" basic model is able to be - retrofitted at time with the „eco" clamping bracket.

#### Can series

The wall brackets are made of stainless steel or powder coated steel. A wall bracket is available for each of the two can sizes. For the can XXL there are 2 variantes available.

#### Exclusive wall bracket

The wall bracket is made of stainless steel. The bucket is compatible with all bucket sizes by adjusting the adjustable mechanism.

## Areas of application

The available different wall bracket variants allow for a site- matched bracket.

## Delivery forms

Wall bracket	Variants	Article no.	Dispenser system
eco	basic	000126	bucket „mini“
		000134	bucket „midi“
		000156	bucket „maxi“ & bucket „XXL“
	safe	000127	bucket „mini“
		000135	bucket „midi“
		000157	bucket „maxi“ & bucket „XXL“
can	clamping bracket		000159
	mini	000317-2	can „mini“
	XXL	000317-3	can „XXL“
		000317-4	can „XXL“
exklusive	universal	000200	all buckets

## Eco series



## Can series



## Exclusive wall bracket



# MediWash

## Hand soap

### Preparation

MediWash is a liquid cream soap for gentle hand washing. Due to the mild substances and skin-friendly pH-value, MediWash is suitable for hand and skin cleansing. The hand soap contains high percentages of skin-nourishing re-fattening agents, mildly readily biodegradable substances and preservatives permitted for cosmetics.

### Application

Lightly moisten your hands and rub well with MediWash. Rinse hands thoroughly.

## Product status / listings

The product complies with the current version of the cosmetics regulation and the legal requirements for biodegradability. The used surfactants fulfil the conditions of biodegradability as specified in Regulation (EC) No. 648/2004 on detergents. MediWash is registered with the EU Cosmetic Products Notification Portal.

## Delivery forms

Name	Containers	Article no.	VPE
MediWash apricot	500 ml	40043	12 bottles
MediWash apricot	10 Litre	40034	1 canister
MediWash apple/kiwi	500 ml	40044	12 bottles
MediWash apple/kiwi	10 Litre	40045	1 canister



# Index

<b>A</b>	
alcohol-free disinfection.....	12/13
alcohol-containing disinfection.....	10/11
<b>B</b>	
buckets.....	20/21
<b>C</b>	
cans.....	20/21
clean & care.....	30-33
cleaningwipes.....	28/29
<b>D</b>	
drywipes.....	6-9
disinfection.....	10-13
<b>F</b>	
FloorWipes.....	26/27
<b>H</b>	
hand cleaning wipes.....	28/29
<b>I</b>	
industry rolls.....	30-33
<b>M</b>	
MediSept.....	10-13
MediWash.....	36/37
MediWipes.....	6-9
<b>N</b>	
nonwoven wipes.....	6-9
<b>P</b>	
PowerWipes.....	28/29
<b>Q</b>	
Quick & Easy.....	16-19
<b>S</b>	
Sachets.....	14/15
saturated disinfection rolls.....	14-19
saturated hand cleaning rolls.....	28/29
soap.....	36/37
SmartBag.....	24/25
SmartPoint.....	22/23
<b>W</b>	
wall brackets.....	34/35





WHR GmbH

Am Galgenturm 2  
97638 Mellrichstadt

Telefon: 0049 (0) 97 76. 70 49 1- 0  
Telefax: 0049 (0) 97 76. 70 49 1- 29

Email: [info@whrgmbh.de](mailto:info@whrgmbh.de)  
Web: [www.whrgmbh.de](http://www.whrgmbh.de)